

ERNST HINRICHS Dental

HinriMill

e4
e5

information



HinriMill .e series

Discover the benefits of CAD/CAM simplicity.

Discover the **benefits**



of CAD/CAM **simplicity**



HinriMill .e4




HinriMill e4 - Perfection for the practice lab.

The e4 is the perfect choice for practice laboratories and sets new standards in efficiency and user-friendliness. The machine is specially designed for the fabrication of chairside ceramic restorations and allows you to quickly and reliably restore your patients in just one session.

The e4 is the ideal entry into the practice laboratory and is particularly suitable for those who want to enter the world of CAD/CAM technology. With the HinriMill e4, you save time and money while improving the quality of your restorations. The machine is extremely easy to operate and the results are always of the highest quality.

Thanks to its compact design and quiet operation, the HinriMill e4 fits into any dental practice. With the HinriMill e4, you are equipped for the future and can offer your patients perfect care every time. Let its performance and precision convince you.


- 4 Axes
- 1 Block
- 7 Tools
- Wet + Dry



The Hinrimill e4 is the ideal solution for fast chairside applications.



*Stefan Tiede,
Dental technician at ERNST HIRNCHS Dental*



HinriMill .e5



Discover the future of milling technology with the HinriMill e5! The e5 is the perfect choice for those who want to achieve high quality and precise results. The machine is easy to operate and enables use in a wide variety of applications with ease.

Thanks to its freedom from compressed air, the e5 offers maximum flexibility in the choice of installation location while reducing your operating costs. The open system makes it easy to get started with the digital production of dental restorations and enables seamless integration into your workflows.

- 5 Axes
- 1 Disc
- 6 Blocks
- 17 Tools



With the automatic tool changer, you can accommodate up to 16 standard tools plus an AIRTOOL for fast and precise machining of blanks and blocks. DENTALCAM CAM software with DIRECT-MILL Technology is included for easy and intuitive operation without additional licensing fees.


The HinriMill e5 is not only compact and lightweight, but also robust and made of high-quality industrial components. As a result, it offers you first-class results with a maximum of efficiency.



With the HinriMill e5 we can mill quickly and precisely without compressed air.



*Andreas Gartner,
Owner of the dental laboratory Gartner in Goslar*



Learn more about HinriMill.e series

Airtool

no compressed air

Discover AIRTOOL - a breakthrough technology that gives the e4 and e5 incredible freedom and flexibility when processing workpieces.

Thanks to the AIRTOOL, our machines no longer require an external compressed air connection or built-in compressor. Instead, the AIRTOOL's turbine blades generate a powerful airflow that keeps the workpiece free of chips and dust during the milling process.

This not only provides a better working environment, but also increases machining efficiency. Join the progress!



innovation

ERNST HINRICHS Dental

Your competent dental partner.



All benefits of .e series



Innovative



Reliable



Fast & Precise





Independent



Economical

ERNST HINRICHS Dental GmbH
Borsigstraße 1, 38644 Goslar

info@hinrichs-dental.de
www.hinrichs-dental.com

 ERNSTHINRICHSdental
 hinrichs_dental